



UltraAmp™

Capture Sequence Cap35 / Oyster®-650 Labeled (50)
3DNA® Dendrimer Signal Amplification Reagent

PRODUCT INFORMATION

Catalog Number:	AC0560
Quantity of 3DNA dendrimer per vial:	1 microgram
Targeting moiety bound to 3DNA dendrimer:	Cap35 Oligonucleotide
Capture sequence on dendrimer:	5' -TCC AAT AGA ATC ACA TCG CTT ACA AGG CAA T -3'
Label moiety bound to 3DNA dendrimer:	Oyster-650 (~50 per dendrimer)
Dendrimer concentration:	10 nanograms of DNA per microliter
Storage conditions:	2-8° Celsius
Shelf life:	Six (6) months from the date of manufacture
3DNA Dendrimer Molecular weight:	1.4 x 10 ⁶ Daltons

DESCRIPTION

A DNA dendrimer reagent containing Cap35 capture sequence oligonucleotide (2-4 per dendrimer molecule) and labeled with ~ 50 Oyster-650 moieties.

STORAGE BUFFER

This UltraAmp reagent is prepared and stored in a buffered solution containing highly pure ethylene glycol, thimerosal (anti-microbial) and other proprietary stabilizer agents.

USAGE

This reagent is effective for the detection of nucleic acid and protein targets containing the complementary nucleic acid sequence to the Cap35 capture sequence. Useful in a variety of assay formats, including DNA and protein microarrays, ELISA and bead based flow assay formats. A quantity of 10 – 50 nanograms of dendrimer reagent per reaction is recommended; this may vary dependent upon the assay type and conditions. We recommend titrating the UltraAmp reagent for each specific assay format to determine the most appropriate quantity required per reaction. UltraAmp reagents are compatible with most assay buffers. Volume excluders, such as dextran sulfate, may be used at concentrations of up to 3% to enhance the performance of these reagents.

APPLICATION NOTES

Application notes detailing the specific use of this product in various assay formats are available at www.genisphere.com/signal_amp_app_notes.html

QUALIFICATION

This reagent has been qualified in a standard DNA microarray assay by binding of the UltraAmp reagent to features containing control DNA probes. This assay demonstrates the ability of the UltraAmp reagent to bind appropriately to its target and provide quantifiable signal. The graph indicates the pixel intensity signals achieved by this reagent and other reagents containing varying quantities of fluorescent labels. This quality control test may not represent the performance of this reagent in this and other types of binding assays.

For research use only. Not for use in diagnostic or therapeutic procedures. Not for resale.

3DNA, Genisphere, and UltraAmp are trademarks of Datascope Corp.

Oyster is a trademark of Denovo Biolabels GmbH.

Protected by U.S. Patents 5,175,270; 5,484,904; 5,487,973, 6,046,038, 6,072,043, 6,110,687 and other patents pending in the U.S. and other countries.

Use of this product may require additional licenses depending on the application.
Genisphere Technical Support Tel: 877.888.3DNA Web: www.genisphere.com
AC0560-V1.01

