



UltraAmp™

Anti-Mouse IgG (Fc) Biotin Labeled (50)
3DNA® Dendrimer Signal Amplification Reagent

PRODUCT INFORMATION

Catalog Number:	AM0570
Quantity of 3DNA dendrimer per vial:	3 micrograms
Antibody bound to 3DNA dendrimer:	Anti-Mouse IgG (Fc) (polyclonal)
Label moiety bound to 3DNA dendrimer:	Biotin (~50 per dendrimer)
Dendrimer concentration:	10 nanograms of DNA per microliter
Storage conditions:	2-8° Celsius
Shelf life:	Six (6) months from the date of manufacture
3DNA Dendrimer Molecular weight:	1.4 x 10 ⁶ Daltons

DESCRIPTION

A DNA dendrimer reagent containing anti-Mouse IgG (Fc) antibody (~ 2-4 antibodies per dendrimer molecule) and labeled with ~ 50 biotin moieties.

STORAGE BUFFER

This UltraAmp reagent is prepared and stored in a buffered solution containing highly pure ethylene glycol, thimerosal (anti-microbial) and other proprietary stabilizer agents.

USAGE

This reagent is effective for the detection of mouse IgG in a variety of assay formats, including ELISA and bead based flow assay formats. A quantity of 25 - 50 nanograms of dendrimer reagent per reaction is recommended; this may vary dependent upon the assay type and conditions. We recommend titrating the UltraAmp reagent for each specific assay format to determine the most appropriate quantity required per reaction. UltraAmp reagents are compatible with most assay buffers. Volume excluders, such as dextran sulfate, may be used at concentrations of up to 3% to enhance the performance of these reagents.

APPLICATION NOTES

Application notes detailing the specific use of this product in various assay formats are available at www.genisphere.com/signal_amp_app_notes.html

QUALIFICATION

This reagent has been qualified in a standard ELISA assay where varying quantities of immobilized mouse IgG were bound by this UltraAmp reagent followed by streptavidin-HRP detection. This assay demonstrates the ability of the UltraAmp reagent to bind appropriately to its target and provide quantifiable signal. The graph indicates the standard curve derived from the quality control testing data and does not represent the performance of this reagent in a typical sandwich immunoassay.

For research use only. Not for use in diagnostic or therapeutic procedures. Not for resale.

3DNA and Genisphere are registered trademarks of Datascope Corp.

UltraAmp is a trademark of Datascope Corp.

Protected by U.S. Patents 5,175,270; 5,484,904; 5,487,973, 6,046,038, 6,072,043, 6,110,687 and other patents pending in the U.S. and other countries.

Use of this product may require additional licenses depending on the application.
Genisphere Technical Support Tel: 877.888.3DNA Web: www.genisphere.com
AM0570-V1.01 © 2000-2005 Genisphere Inc., all rights reserved.

