



## Two-Round Linear Amplification Kit

You have 10ng of RNA that looks like this

And you want to compare it to RNA that looks like this



RNA isolated from 5 year old Formalin-Fixed Paraffin Embedded (FFPE) colon cancer tissue sample.



RNA isolated from a biologically similar fresh frozen colon cancer sample.

RampUP amplified senseRNA prepared from the above samples demonstrated a Pearson correlation of 0.902 when compared on Affymetrix U133A GeneChips.

### RampUP's TRPP Technology Provides:

- Simple protocol
- Amplification of any RNA
- Artifact-free amplification
- Sense-strand amplified RNA
- Up to 10<sup>6</sup>-fold amplification
- Accurate and reproducible results

**RampUP produces senseRNA that is adaptable to any RNA expression assay!**

### Ordering Information:

Product	Number of Assays	Catalog Number
RampUP	10	RAMPUP10
RampUP Plus™	10	RPLUS10
Plus	10	P10PLUS
cDNA Synthesis	10	CDNAMOD

**Genisphere®**

SIGNAL + SAMPLE AMPLIFICATION PRODUCTS

### To Order:

Genisphere Inc.  
 2801 Sterling Drive  
 Hatfield, PA 19440  
 Tel: 877.888.3DNA  
 215.996.3002  
 Fax: 877.FAX.3DNA  
 215.996.3070  
[www.genisphere.com](http://www.genisphere.com)